

# Colors of Manganese Ions

## Purpose

To illustrate the various oxidation states of an element.

## Materials

- Potassium permanganate,  $\text{KMnO}_4$  (0.1 M)
- Sodium hydroxide,  $\text{NaOH}$  (2.0 M)
- Sulfuric acid,  $\text{H}_2\text{SO}_4$  (3.0 M)
- Sodium bisulfite,  $\text{NaHSO}_3$  (0.1 M)
- EDTA (1.0 M)
- Manganese(II) sulfate,  $\text{MnSO}_4$  (powder)
- (5) 100-mL beakers
- 50 mL graduated cylinder

## Safety

- Read the SDS sheets for all chemicals before using them.
- Wear safety glasses, gloves, and lab coat.
- Concentrated acids and bases are used.
- Permanganate solution will stain.

## Procedure

1. Pour 50 mL of 0.1 M potassium permanganate in five 100-mL beakers (labeled 1 to 5).
2. To beaker number 1 add 15 mL of 3.0 M sulfuric acid and then, while stirring slowly, add 0.1 M sodium bisulfite until a color change takes place.
3. To beaker number 2 add 20 mL of 2.0 M sodium hydroxide and then, while stirring slowly, add 0.1 M sodium bisulfite until a color change takes place.
4. To beaker number 3 slowly add 0.1 M sodium bisulfite while stirring until a color change takes place.
5. To beaker number 4 add 5 mL of 1.0 M EDTA and a pinch of solid manganese(II) sulfate; then stir.

## Follow-up Teaching Notes

- The chart below outlines color corresponding to the various oxidation states of manganese

Beaker	Oxidation State of Mn	Color
1	+2	colorless
2	+6	green
3	+4	brown
4	+3	violet/rose
5	+7	purple

## Concepts

- Oxidation states, balancing redox equations.

## Extension

- Students can be asked to write balanced equations for beakers 1, 2, and 3 if they are told the oxidation half-reaction is  $\text{HSO}_3^{-1} \rightarrow \text{SO}_4^{-2}$
- The reduction half-reaction is  $\text{MnO}_4^{-1} \rightarrow \text{Mn}^?$  where the ? refers to the oxidation state of manganese in the appropriate beaker.

## Disposal/Clean-up

- Remaining contents can be placed in the science department's heavy metal waste container for proper disposal.